

# Product datasheet

## G2 Duoplus

*This treatment couch offers the ideal solution for treatment and massage of patients in a lying position. The head rest can be adjusted upward and downward extremely easily with a gas spring. The G2 Duoplus treatment and examination couch can be equipped with armrests for optimal relaxation in the prone position.*

### Standard features

- Two-section couch with adjustable head section
- Head section supported by gas spring
- Electric height adjustment with foot switch
- Cushions with 40 mm padding
- 15 cm clearance underneath frame for use with patient hoist (not compatible with 150 007)
- Non slip feet
- Cushions in a choice of 11 colours

### Ordering data

160 020 G2 Duoplus

### Optional

- 150 029 G2 Hand control
- 150 007 G2 Peripheral switch
- 150 001 G2 Hydraulic
- 150 062 G2 Wheel mechanism
- 150 002 G2 Arm supports\*
- 150 003 G2 Cushion with face opening
- 150 004 G2 Cushions in non-standard colour

### Accessories

- 150 106 G2 Paper roll holder 50 cm
- 150 117 G2 Paper roll holder 60 cm
- 337 018 Paper rolls 50 cm (per 12 pieces)
- 332 915 Foot end protector for cushion
- xxx xxx Terry-cloth cover, stretch material (available in 9 colours 'Style')

\* Automatically includes 150 003

### Specifications

Length x width:	± 205 x 67cm
Height adjustment:	± 47 – 100 cm
Weight:	± 80 kg
Max. lifting capacity:	200 kg
Electric motor:	6.000 N
Hydraulic pump:	10.000 N
Power supply:	230 VAC, 50-60 Hz, 1,0 A 120 VAC, 60 Hz, 1,7 A
Insulation class:	II
Moister resistance:	IP 21
Classification:	I
Conformity:	IEC/EN 60601-1 IEC/EN 60601-2 MDD 93/42/EEC
<b>CE</b>	
Guarantee:	3 years (upholstery 1 year)




## G2 Duoplus



G2 Duoplus shown with options:

- Hydraulic height adjustment
- Wheel mechanism
- Arm supports

### Functionality/dimensions



### Cushion colours

#### Style



#### Fashion

