

# Product datasheet

## Osteoflex Advanced

Particularly suitable for those who concentrate on treatment of the spinal column. Thanks to the unique fold-down side supports, you can stand closer to the patient. This enables you to exert more force with less strain on your own shoulders and back. The arm supports provide a very comfortable ventral position, with the best possible relaxation of the shoulder and neck muscles. Perfectly suited for osteopathic manipulations.

### Standard features

- Three-section couch with adjustable head-middle- and leg section
- All sections gas spring supported
- Arm supports in the head section
- Unique fold-down supports in the middle section
- Electric or hydraulic height adjustment
- Gymna Stability Profile (GSP)
- Cushions with 40 mm padding
- Breathing hole in head section, incl. cushion
- Central wheel-lifting mechanism with 4 casters
- Non slip feet
- Cushions in a choice of 11 colours
- Prepared for hip traction support
- Prepared for shoulder traction supports

### Ordering data

- 125.061 Osteoflex Advanced peripheral switch  
 125 062 Osteoflex Advanced foot switch  
 125 063 Osteoflex Advanced hydraulic

### Specifications

Length x width:	± 202 x 67cm
Height adjustment:	± 45 – 95 cm
Weight:	± 105 kg
Max. lifting capacity:	200 kg
Electric motor:	6.000 N
Hydraulic pump:	10.000 N
Power supply:	100 – 240 VAC, 50-60 Hz
Current consumption:	1,5 A
Insulation class:	II
Moister resistance:	IP 21
Classification:	I
Conformity:	IEC/EN 60601-1 IEC/EN 60601-2 CAN/CSA C22.2 No. 601.1 MDD 93/42/EEC
Guarantee:	3 years (upholstery 1 year)



### Optional

- 125 602 Electrical adjustable middle section  
 125 610 Motorized head section (not in combination with 125 602)  
 121 405 Set brackets for fastening straps  
 121 402 Cushions 60 cm wide (Advanced)

### Accessories

- 122 405 Paper roll holder, for undercarriage  
 337 018 Paper rolls (per 12 pieces)  
 121 409 Extra footswitch for height adjustment  
 125 652 Hand switch for height adjustment  
 107 819 Traction support for hips or shoulders  
 121 401 Gynaecological leg support  
 333 245 Traction device AutoTrac 460  
 121 434 Attachment frame for AutoTrac 450/460  
 332 915 Foot end protector for cushion  
 xxx xxx Terry-cloth cover, stretch material (available in 9 colours 'Style')





Gymna Stability Profile

**Osteoflex Advanced**



**Functionality/dimensions**



**Cushion colours**

**Style**



**Fashion**

