

Product datasheet

Trioflex Advanced

The standard treatment couch for physiotherapy! Adjustable headrest with arm supports, middle-and leg section. Extremely versatile and suited for treatment in sitting, lying, lumbar flexion and Trendelenburg position. Guarantees years of problem-free support for your manipulations.

The arm supports provide a very comfortable ventral position, with the best possible relaxation of the shoulder and neck muscles.

Standard features

- Three-section couch with adjustable head-, middle-, and leg section
- All sections gas spring supported
- Arm supports in the head section
- Electric or hydraulic height adjustment
- Gymna Stability Profile (GSP)
- Cushions with 40 mm padding
- Breathing hole in head section, incl. cushion
- Central wheel-lifting mechanism with 4 casters
- Non slip feet
- Cushions in a choice of 11 colours
- Prepared for hip traction support
- Prepared for shoulder traction supports

Ordering data

125 051 Trioflex Advanced peripheral switch

125 052 Trioflex Advanced foot switch125 053 Trioflex Advanced hydraulic

Specifications

Hydraulic pump:

Length x width: $\pm 202 \times 67 \text{cm}$ Height adjustment: $\pm 45 - 95 \text{ cm}$ Weight: $\pm 100 \text{ kg}$ Max. lifting capacity: $\pm 200 \text{ kg}$ Electric motor: $\pm 6.000 \text{ N}$

Power supply: 100 – 240 VAC, 50-60 Hz

10.000 N

Current consumption: 1,5 A
Insulation class: II
Moister resistance: IP 21
Classification: I

Conformity: IEC/EN 60601-1

IEC/EN 60601-2
CAN/CSA C22.2 No. 601.1

MDD 93/42/EEC

Guarantee: 3 years (upholstery 1 year)

Optional

125 602 Electrical adjustable middle section
125 610 Motorized head section (not in combination with 125 602)
121 404 Set brackets for fastening straps
121 402 Cushions 60 cm wide (Advanced)
121 403 Cushions 80 cm wide (Advanced)

Accessories

122 405 Paper roll holder, for undercarriage 337 018 Paper rolls (per 12 pieces) 121 409 Extra footswitch for height adjustment 125 652 Hand switch for height adjustment 107 819 Traction support for hips or shoulders 121 401 Gynaecological leg support 333 245 Traction device AutoTrac 460 121 434 Attachment frame for AutoTrac 450/460 332 915 Foot end protector for cushion Terry-cloth cover, stretch material XXX XXX (available in 9 colours 'Style')



Gymna Stability Profile

gymna

Trioflex Advanced



Functionality/dimensions -72° 0-75° 46 54 102 Option 60 cm Option 80 cm

Cushion colours

